

TOP GLOVE SDN BHD

INSTRUCTION FOR USE

EA101 - LATEX EXAMINATION POWDERED GLOVES

Sphere of Use: Latex Examination Powdered gloves are worn to protect the hand of users from chemical and also to prevent transmission of diseases.

Instructions : **Caution :-** User should consider circumstance of use in after donning. Powder can be removed by wiping gloves thoroughly with a sterile wet sponge, sterile wet towel, or other effective methods. "Safe use of this gloves by or on latex sensitized individuals has not been established".


Warning :- Isolated cases of allergic reactions to rubber latex or powder have been reported. If you experience a reaction to this product, discontinue use immediately and consult your physician. The packing of this product contains Natural Rubber Latex which may cause allergic reaction in some individuals.

"Other components used in making gloves may cause allergic reactions in some users"


:Storage :- Store in a cool, dry place. Opened box should be shielded from exposure to direct sun or fluorescent lighting.

Shelf Life : 5 years upon the manufacturing date.

Additional Info : Sizes available : XS - XL (EN420:2003+A1:2009)




EU Type Examination and ongoing Conformity by Notified Body:- **CE2777**
SATRA Technology Europe Ltd
Bracetown Business Park,
Clonee, D15 YN2P, Ireland.



Hand Size	6 (XS)	7 (S)	8 (M)	9 (L)	10 (XL)
Min Length	220mm	230mm	240mm	250mm	260mm

Resistance against Bacteria and Fungi - **PASS**
Resistance against Virus - **PASS**

EN ISO 374-1:2016 Classification Of Permeation Performance Level	
Measured Breakthrough time (min)	Permeation Performance Level.
> 10	1
> 30	2
> 60	3
> 120	4
>240	5
>480	6



Code letter	Chemical Permeation (EN ISO 374-1:2016)	Level	Mean Degradation (%) (EN374-4:2013)	Degradation levels indicate the change in Puncture Resistance of the glove after exposure to the tested chemical.
K	40% Sodium Hydroxide	2	13.8	
M	65% Nitric Acid	1	46.8	
T	37% Formaldehyde	2	-5.1	
P	30% Hydrogen Peroxide	1	-16.1	
S	40% Hydrofluoric Acid	2	N/A	

EA101 - STATEMENT AND CAUTION

- 1) This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
- 2) The penetration resistance and chemical resistance have **been assessed under laboratory conditions and relates only to the tested specimen** taken from the palm only (except in cases where the glove is equal to or over 400mm – where the cuff is tested also) and **relates only to the chemical tested**. It can be different if the chemical is used in a mixture.
- 3) It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- 4) When used, protective gloves may provide less resistance to the dangerous chemical due to changes in the physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
- 5) Before usage, inspect the gloves for any defects or imperfections.